

Product Datasheet

Mouse Arginase-2 ELISA Kit (orb1210845)

Description	Mouse Arginase-2, ARG2 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	ARG2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	38,878
Uniprot ID	O08691
Solubility (25°C)	24ng/mL
Sensitivity	0.026ng/mL
Entrez	11847
Expiration Date	Please enquire.