

Product Datasheet

Mouse Anti Beta-1 Adrenergic Receptor ELISA Kit (orb1210848)

Description	Mouse Anti Beta-1 Adrenergic Receptor, ADRB1-AB ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	beta1AR antibody
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	50,494
Uniprot ID	P34971
Solubility (25°C)	24ng/mL
Sensitivity	0.023ng/mL
Entrez	11554
Expiration Date	Please enquire.