

Product Datasheet

Rat Atp Translocase 1 ELISA Kit (orb1210851)

Description	Rat Atp Translocase 1, SLC25A4 ELISA Kit
Reactivity	Rat
Range	10-2000ng/L
Tested Applications	ELISA
Target	SLC25A4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	32,989
Uniprot ID	Q05962
Solubility (25°C)	2400ng/L
Sensitivity	4.23ng/L
Entrez	85333
Expiration Date	Please enquire.