

Product Datasheet

Mouse Adenosine Receptor A2A ELISA Kit (orb1210865)

Description	Mouse Adenosine Receptor A2A, ADORA2A ELISA Kit
Reactivity	Mouse
Range	20-4500ng/L
Tested Applications	ELISA
Target	ADORA2A
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	44,971
Uniprot ID	Q60613
Solubility (25°C)	4800ng/L
Sensitivity	9.41ng/L
Entrez	11540
Expiration Date	Please enquire.