

Product Datasheet

Mouse Beta-defensin 2 ELISA Kit (orb1210912)

Description	Mouse Beta-defensin 2, DEFB2 ELISA Kit
Reactivity	Mouse
Range	5-1000ng/L
Tested Applications	ELISA
Target	DEFB2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,820
Uniprot ID	P82020
Solubility (25°C)	1280ng/L
Sensitivity	2.52ng/L
Entrez	13215
Expiration Date	Please enquire.