

Product Datasheet

Mouse Transcription Factor A ELISA Kit (orb1210931)

Description	Mouse Transcription Factor A, TFAM ELISA Kit
Reactivity	Mouse
Range	5-1500ng/L
Tested Applications	ELISA
Target	TFAM
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,988
Uniprot ID	P40630
Solubility (25°C)	1600ng/L
Sensitivity	2.41ng/L
Entrez	21780
Expiration Date	Please enquire.