

Product Datasheet

Mouse Dihydropyrimidinase Like Protein 2 ELISA Kit (orb1210932)

Description	Mouse Dihydropyrimidinase Like Protein 2, DPYSL2 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	DPYSL2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	62,278
Uniprot ID	008553
Solubility (25°C)	48ng/ml
Sensitivity	0.04ng/ml
Entrez	12934
Expiration Date	Please enquire.