

Product Datasheet

Mouse Adenosine Deaminase ELISA Kit (orb1210939)

Description	Mouse Adenosine Deaminase, ADA ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	ADA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,992
Uniprot ID	P03958
Solubility (25°C)	40ng/ml
Sensitivity	0.023ng/ml
Entrez	11486
Expiration Date	Please enquire.