


Product Datasheet

Mouse Mast Cell Protease 1 ELISA Kit (orb1210942)

Description	Mouse Mast Cell Protease 1, MCPT1 ELISA Kit
Reactivity	Mouse
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	mcpt1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,013
Uniprot ID	P11034
Solubility (25°C)	16ng/ml
Sensitivity	0.01ng/ml
Entrez	17224
Expiration Date	Please enquire.



Standard Curve of Mouse Mast Cell Protease 1 ELISA Kit