

Product Datasheet

Mouse A Disintegrin and Metalloproteinase with Thrombospondin Motifs 8 ELISA Kit (orb1210944)

Description	Mouse A Disintegrin and Metalloproteinase with Thrombospondin Motifs 8, ADAMTS8 ELISA Kit
Reactivity	Mouse
Range	6-380ng/ml
Tested Applications	ELISA
Target	ADAMTS8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	98,880
Uniprot ID	P57110
Solubility (25°C)	400ng/ml
Sensitivity	4.12ng/ml
Entrez	30806
Expiration Date	Please enquire.