

Product Datasheet

Human G Protein-regulated Inducer of Neurite Outgrowth 1 ELISA Kit (orb1210946)

Description	Human G Protein-regulated Inducer of Neurite Outgrowth 1, GPRIN1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	GPRIN1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	102,399
Uniprot ID	Q7Z2K8
Solubility (25°C)	2400ng/L
Sensitivity	5.03ng/L
Entrez	114787
Expiration Date	Please enquire.