

Product Datasheet

Mouse S100 Calcium Binding Protein A7 ELISA Kit (orb1210985)

Description	Mouse S100 Calcium Binding Protein A7, S100A7 ELISA Kit
Reactivity	Mouse
Range	10-4000ng/L
Tested Applications	ELISA
Target	S100A7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	12,870
Uniprot ID	Q6S5I3
Solubility (25°C)	4800ng/L
Sensitivity	5.42ng/L
Entrez	381493
Expiration Date	Please enquire.