

Product Datasheet

Mouse Adropin ELISA Kit (orb1211005)

Description	Mouse Adropin, AD ELISA Kit
Reactivity	Mouse
Range	50-10000ng/L
Tested Applications	ELISA
Target	AD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	7,927
Uniprot ID	Q8K1D8
Solubility (25°C)	12800ng/L
Sensitivity	24.95ng/L
Entrez	69638
Expiration Date	Please enquire.