

Product Datasheet

Rat Nadphoxidase 4 ELISA Kit (orb1211012)

Description	Rat Nadphoxidase 4, NOX4 ELISA Kit
Reactivity	Rat
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	NOX4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	66,469
Uniprot ID	Q924V1
Solubility (25°C)	64ng/ml
Sensitivity	0.11ng/ml
Entrez	85431
Expiration Date	Please enquire.