

Product Datasheet

Human A Disintegrin and Metalloproteinase with Thrombospondin Motifs 7 ELISA Kit (orb1211052)

Description	Human A Disintegrin and Metalloproteinase with Thrombospondin Motifs 7, ADAMTS7 ELISA Kit
Reactivity	Human
Range	1-400ng/ml
Tested Applications	ELISA
Target	ADAMTS7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	184,095
Uniprot ID	Q9UKP4
Solubility (25°C)	480ng/ml
Sensitivity	0.65ng/ml
Entrez	11173
Expiration Date	Please enquire.