

Product Datasheet

Rat Pro-Opiomelanocortin ELISA Kit (orb1211055)

Description	Rat Pro-Opiomelanocortin, POMC ELISA Kit
Reactivity	Rat
Range	3-900ng/L
Tested Applications	ELISA
Target	POMC
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	26,540
Uniprot ID	P01194
Solubility (25°C)	960ng/L
Sensitivity	1.61ng/L
Entrez	24664
Expiration Date	Please enquire.