

Product Datasheet

Rat Glucose Transporter 2 ELISA Kit (orb1211074)

Description	Rat Glucose Transporter 2, GLUT2 ELISA Kit
Reactivity	Rat
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	GLUT2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,086
Uniprot ID	P12336
Solubility (25°C)	80ng/ml
Sensitivity	0.12ng/ml
Entrez	25351
Expiration Date	Please enquire.