

Product Datasheet

Mouse Dystrophin ELISA Kit (orb1211083)

Description	Mouse Dystrophin, DMD ELISA Kit
Reactivity	Mouse
Range	0.2-60ng/ml
Tested Applications	ELISA
Target	DMD
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	425,832
Uniprot ID	P11531
Solubility (25°C)	64ng/ml
Sensitivity	0.097ng/ml
Entrez	13405
Expiration Date	Please enquire.