

Product Datasheet

Rat Hyaluronan Mediated Motility Receptor ELISA Kit (orb1211085)

Description	Rat Hyaluronan Mediated Motility Receptor, HMMR ELISA Kit
Reactivity	Rat
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	HMMR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	57,858
Uniprot ID	P97779
Solubility (25°C)	24ng/ml
Sensitivity	0.016ng/ml
Entrez	25460
Expiration Date	Please enquire.