

Product Datasheet

Human Phospholipid Scramblase 1 ELISA Kit (orb1211088)

Description	Human Phospholipid Scramblase 1, PLSCR1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	PLSCR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	35,049
Uniprot ID	015162
Solubility (25°C)	3200ng/L
Sensitivity	8.05ng/L
Entrez	5359
Expiration Date	Please enquire.