

Product Datasheet

Rat Cytochrome P450 19A1 ELISA Kit (orb1211094)

Description	Rat Cytochrome P450 19A1, CYP19A1 ELISA Kit
Reactivity	Rat
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	CYP19A1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	58,412
Uniprot ID	P22443
Solubility (25°C)	32ng/ml
Sensitivity	0.03ng/ml
Entrez	25147
Expiration Date	Please enquire.