

Product Datasheet

Human L-dopachrome Tautomerase ELISA Kit (orb1211099)

Description	Human L-dopachrome Tautomerase, DCT ELISA Kit
Reactivity	Human
Range	0.2-70ng/ml
Tested Applications	ELISA
Target	DCT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	59,145
Uniprot ID	P40126
Solubility (25°C)	80ng/ml
Sensitivity	0.12ng/ml
Entrez	1638
Expiration Date	Please enquire.