

Product Datasheet

Rat A Disintegrin and Metalloproteinase with Thrombospondin Motifs 1 ELISA Kit (orb1211106)

Description	Rat A Disintegrin and Metalloproteinase with Thrombospondin Motifs 1, ADAMTS1 ELISA Kit
Reactivity	Rat
Range	0.05-15ng/ml
Tested Applications	ELISA
Target	ADAMTS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	105,706
Uniprot ID	Q9WUQ1
Solubility (25°C)	16ng/ml
Sensitivity	0.027ng/ml
Entrez	79252
Expiration Date	Please enquire.