

Product Datasheet

Mouse Caspase 12 ELISA Kit (orb1211121)

Description	Mouse Caspase 12, CASP12 ELISA Kit
Reactivity	Mouse
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	CASP12
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	47,854
Uniprot ID	O08736
Solubility (25°C)	24ng/ml
Sensitivity	0.026ng/ml
Entrez	12364
Expiration Date	Please enquire.