

Product Datasheet

Human Guanine Nucleotide Exchange Factor Vav3 ELISA Kit (orb1211143)

Description	Human Guanine Nucleotide Exchange Factor Vav3, VAV3 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	VAV3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	97,776
Uniprot ID	Q9UKW4
Solubility (25°C)	24ng/ml
Sensitivity	0.02ng/ml
Entrez	10451
Expiration Date	Please enquire.