

Product Datasheet

Human Latent Transforming Growth Factor Beta Binding Protein 2 ELISA Kit (orb1211198)

Description	Human Latent Transforming Growth Factor Beta Binding Protein 2, LTBP2 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	LTBP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	195,052
Uniprot ID	Q14767
Solubility (25°C)	24ng/ml
Sensitivity	0.023ng/ml
Entrez	4053
Expiration Date	Please enquire.