

Product Datasheet

Bovine Ubiquitin-like Protein Isg15 ELISA Kit (orb1211219)

Description	Bovine Ubiquitin-like Protein Isg15, ISG15 ELISA Kit
Reactivity	Bovine
Range	5-2000ng/L
Tested Applications	ELISA
Target	ISG15
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,310
Uniprot ID	002741
Solubility (25°C)	2400ng/L
Sensitivity	2.45ng/L
Entrez	281871
Expiration Date	Please enquire.