

Product Datasheet

Human Polymeric Immunoglobulin Receptor Secretory Component ELISA Kit (orb1211241)

Description	Human Polymeric Immunoglobulin Receptor Secretory Component, PIGR ELISA Kit
Reactivity	Human
Range	1-300ng/ml
Tested Applications	ELISA
Target	PIGR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	83,284
Uniprot ID	P01833
Solubility (25°C)	320ng/ml
Sensitivity	0.51ng/ml
Entrez	5284
Expiration Date	Please enquire.