

Product Datasheet

Human Vitronectin ELISA Kit (orb1211257)

Description	Human Vitronectin, VTN ELISA Kit
Reactivity	Human
Range	10-3000ng/L
Tested Applications	ELISA
Target	VTN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	54,306
Uniprot ID	P04004
Solubility (25°C)	3200ng/L
Sensitivity	5.12ng/L
Entrez	7448
Expiration Date	Please enquire.