

Product Datasheet

Chicken Ovotransferrin ELISA Kit (orb1211258)

Description	Chicken Ovotransferrin, OTF ELISA Kit
Reactivity	Gallus
Range	0.1-35ug/ml
Tested Applications	ELISA
Target	OTF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	77,777
Uniprot ID	P02789
Solubility (25°C)	40µg/ml
Sensitivity	0.047µg/ml
Entrez	396241
Expiration Date	Please enquire.