

Product Datasheet

Chicken Bone Morphogenetic Protein 2 ELISA Kit (orb1211269)

Description	Chicken Bone Morphogenetic Protein 2, BMP2 ELISA Kit
Reactivity	Gallus
Range	20-4500ng/L
Tested Applications	ELISA
Target	BMP2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,347
Uniprot ID	Q90751
Solubility (25°C)	4800ng/L
Sensitivity	9.78ng/L
Entrez	378779
Expiration Date	Please enquire.