

Product Datasheet

Human Leukocyte Immunoglobulin-like Receptor Subfamily B Member 4 ELISA Kit (orb1211275)

Description	Human Leukocyte Immunoglobulin-like Receptor Subfamily B Member 4, LILRB4 ELISA Kit
Reactivity	Human
Range	6-2000ng/L
Tested Applications	ELISA
Target	LILRB4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	49,356
Uniprot ID	Q8NHJ6
Solubility (25°C)	2400ng/L
Sensitivity	3.22ng/L
Entrez	11006
Expiration Date	Please enquire.