

Product Datasheet

Human C-type Lectin Domain Family 4 Member E ELISA Kit (orb1211280)

Description	Human C-type Lectin Domain Family 4 Member E, CLEC4E ELISA Kit
Reactivity	Human
Range	5-2000ng/L
Tested Applications	ELISA
Target	CLEC4E
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	25,073
Uniprot ID	Q9ULY5
Solubility (25°C)	2400ng/L
Sensitivity	2.33ng/L
Entrez	26253
Expiration Date	Please enquire.