

Product Datasheet

Sheep Insulin-like Growth Factor 2 ELISA Kit (orb1211305)

Description	Sheep Insulin-like Growth Factor 2, IGF2 ELISA Kit
Reactivity	Sheep
Range	6-1800ng/mL
Tested Applications	ELISA
Target	IGF2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,616
Uniprot ID	P10764
Solubility (25°C)	1920ng/mL
Sensitivity	2.68ng/mL
Entrez	443325
Expiration Date	Please enquire.