

Product Datasheet

Rat Interleukin-7 ELISA Kit (orb1211323)

Description	Rat Interleukin-7, IL-7 ELISA Kit
Reactivity	Rat
Range	10-2000ng/L
Tested Applications	ELISA
Target	IL-7
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,709
Uniprot ID	P56478
Solubility (25°C)	2400ng/L
Sensitivity	5.24ng/L
Entrez	25647
Expiration Date	Please enquire.