

Product Datasheet

Sheep Insulin-like Growth Factors 1 ELISA Kit (orb1211343)

Description	Sheep Insulin-like Growth Factors 1, IGF-1 ELISA Kit
Reactivity	Sheep
Range	5-1500ng/ml
Tested Applications	ELISA
Target	IGF-1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,012
Uniprot ID	P10763
Solubility (25°C)	1600ng/ml
Sensitivity	2.39ng/ml
Entrez	443318
Expiration Date	Please enquire.