

Product Datasheet

Human Interleukin 5 ELISA Kit (orb1211346)

Description	Human Interleukin 5, IL-5 ELISA Kit
Reactivity	Human
Range	3-900ng/L
Tested Applications	ELISA
Target	IL-5
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	15,238
Uniprot ID	P05113
Solubility (25°C)	960ng/L
Sensitivity	1.52ng/L
Entrez	3567
Expiration Date	Please enquire.