

Product Datasheet

Horse Brain Derived Neurotrophic Factor ELISA Kit (orb1211347)

Description	Horse Brain Derived Neurotrophic Factor, BDNF ELISA Kit
Reactivity	Equine
Range	5-2000ng/L
Tested Applications	ELISA
Target	BDNF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,722
Uniprot ID	Q0EAB7
Solubility (25°C)	2400ng/L
Sensitivity	2.92ng/L
Entrez	627
Expiration Date	Please enquire.