

Product Datasheet

Cat Trypsin ELISA Kit (orb1211361)

Description	Cat Trypsin, TRP ELISA Kit
Reactivity	Feline
Range	0.5-200ng/ml
Tested Applications	ELISA
Target	Trypsin
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	1,825
Uniprot ID	P81071
Solubility (25°C)	240ng/ml
Sensitivity	0.28ng/ml
Entrez	101093891
Expiration Date	Please enquire.