

Product Datasheet

Sheep Interleukin 6 ELISA Kit (orb1211364)

Description	Sheep Interleukin 6, IL-6 ELISA Kit
Reactivity	Sheep
Range	2-600ng/L
Tested Applications	ELISA
Target	IL-6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,446
Uniprot ID	P29455
Solubility (25°C)	640ng/L
Sensitivity	1.06ng/L
Entrez	443406
Expiration Date	Please enquire.