

Product Datasheet

Sheep Leptin ELISA Kit (orb1211366)

Description	Sheep Leptin, LEP ELISA Kit
Reactivity	Sheep
Range	1-400ng/mL
Tested Applications	ELISA
Target	LEP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	16,054
Uniprot ID	Q28603
Solubility (25°C)	480ng/mL
Sensitivity	0.52ng/mL
Entrez	443534
Expiration Date	Please enquire.