

Product Datasheet

Goat Leptin ELISA Kit (orb1211367)

Description	Goat Leptin, LEP ELISA Kit
Reactivity	Goat
Range	0.5-100ng/ml
Tested Applications	ELISA
Target	LEP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,700
Uniprot ID	Q257X2
Solubility (25°C)	128ng/ml
Sensitivity	0.27ng/ml
Entrez	100860794
Expiration Date	Please enquire.