

Product Datasheet

Cat Interferon Beta ELISA Kit (orb1211377)

Description	Cat Interferon Beta, IFN-beta ELISA Kit
Reactivity	Feline
Range	2-700ng/L
Tested Applications	ELISA
Target	IFN-beta
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,188
Uniprot ID	Q9N2J0
Solubility (25°C)	800ng/L
Sensitivity	1.22ng/L
Entrez	493849
Expiration Date	Please enquire.