

Product Datasheet

Sheep Interleukin 1Beta ELISA Kit (orb1211378)

Description	Sheep Interleukin 1Beta, IL-1B ELISA Kit
Reactivity	Sheep
Range	0.5-150ng/L
Tested Applications	ELISA
Target	IL-1beta
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,717
Uniprot ID	P21621
Solubility (25°C)	160ng/L
Sensitivity	0.24ng/L
Entrez	443539
Expiration Date	Please enquire.