

## Product Datasheet

### Human Insulin Receptor Substrate 2 ELISA Kit (orb1211384)

|                                |                                                        |
|--------------------------------|--------------------------------------------------------|
| <b>Description</b>             | Human Insulin Receptor<br>Substrate 2, IRS-2 ELISA Kit |
| <b>Reactivity</b>              | Human                                                  |
| <b>Range</b>                   | 0.5-200ng/L                                            |
| <b>Tested<br/>Applications</b> | ELISA                                                  |
| <b>Target</b>                  | IRS-2                                                  |
| <b>Storage</b>                 | 2-8°C                                                  |
| <b>Note</b>                    | For research use only                                  |
| <b>Application notes</b>       | Sandwich                                               |
| <b>Sample Types</b>            | Serum, plasma, cell culture<br>supernates              |
| <b>MW</b>                      | 137,334                                                |
| <b>Uniprot ID</b>              | <b>Q9Y4H2</b>                                          |
| <b>Solubility (25°C)</b>       | 240ng/L                                                |
| <b>Sensitivity</b>             | 0.23ng/L                                               |
| <b>Entrez</b>                  | <b>8660</b>                                            |
| <b>Expiration Date</b>         | Please enquire.                                        |