

Product Datasheet

Cat Interleukin 10 ELISA Kit (orb1211390)

Description	Cat Interleukin 10, IL-10 ELISA Kit
Reactivity	Feline
Range	20-4000ng/L
Tested Applications	ELISA
Target	IL-10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,575
Uniprot ID	P55029
Solubility (25°C)	4800ng/L
Sensitivity	10.51ng/L
Entrez	493683
Expiration Date	Please enquire.