

Product Datasheet

Chicken Leptin ELISA Kit (orb1211394)

Description	Chicken Leptin, LEP ELISA Kit
Reactivity	Gallus
Range	5-1500ng/ml
Tested Applications	ELISA
Target	LEP
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	18,183
Uniprot ID	O42164
Solubility (25°C)	1920ng/ml
Sensitivity	2.43ng/ml
Entrez	373955
Expiration Date	Please enquire.