

Product Datasheet

Mouse Angiotensin 2 ELISA Kit (orb1211400)

Description	Mouse Angiotensin 2, ANG-2 ELISA Kit
Reactivity	Mouse
Range	0.2-60ng/L
Tested Applications	ELISA
Target	ANG-2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,990
Uniprot ID	P11859
Solubility (25°C)	64ng/L
Sensitivity	0.124ng/L
Entrez	11606
Expiration Date	Please enquire.