

Product Datasheet

Rat Prostaglandin-H2 D-isomerase ELISA Kit (orb1211408)

Description	Rat Prostaglandin-H2 D-isomerase, PTGDS ELISA Kit
Reactivity	Rat
Range	5-1000ng/mL
Tested Applications	ELISA
Target	PTGDS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,301
Uniprot ID	P22057
Solubility (25°C)	1280ng/mL
Sensitivity	2.14ng/mL
Entrez	25526
Expiration Date	Please enquire.