

Product Datasheet

Rat Muscarinic Acetylcholine Receptor M1 ELISA Kit (orb1211428)

Description	Rat Muscarinic Acetylcholine Receptor M1, CHRM1 ELISA Kit
Reactivity	Rat
Range	10-2000ng/L
Tested Applications	ELISA
Target	CHRM1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	51,369
Uniprot ID	P08482
Solubility (25°C)	2400ng/L
Sensitivity	5.46ng/L
Entrez	25229
Expiration Date	Please enquire.